

OL 770 Multi-channel Spectroradiometer Recommended Computer Configuration

	Optimum Configuration	Minimum Configuration
	For optimum performance, use the following configuration:	The following configuration is adequate for OL 770 operation, but there will be noticeable decreases in response time and system performance:
Computer:	1.5 GHz Pentium III® or Athlon® -based IBM® - compatible personal computer or laptop	Pentium® II- based IBM® -compatible personal computer or laptop, 500 MHz
Memory:	64 Mb RAM or greater	32 Mb RAM
Operating System:	Windows® 2000/XP	SAME
Monitor:	SVGA capable of displaying 1024x768 or greater in 24 bit color	SVGA capable of displaying 800x600 in 16 bit color
Mass Storage:	Hard Drive > 1 Gb	> 20 Mb free space
Mouse:	Microsoft™ compatible	SAME
Printer:	Windows® supported printer	SAME
Interface (standard):	Serial port - RS232, USB	RS232 ~ 2 spectra/sec
	USB port and Windows® 2000 or XP	SAME
Software:	Microsoft™ Excel® & Word®	SAME

* Use of the minimum system configuration will limit transfer rates in USB mode.

